



ROOTT

ROOTT 



Meet intelligence solutions with

Direct scan



Time efficiency

No traditional impression, no need to install scan posts



Cost-efficiency

Less chair-time, no extra parts, scan bodies




Easy & smart

More convenience with less patient operations,
less components to handle

Prosthetic workflow

ROOTT **C** **CS**



ROOTT **C**

Cement & telescopic retained


One-piece implant

Simple solution to bone atrophy


One-piece implant for more comfort and simplicity with a bendable neck for up to 15°. It ensures extreme time and cost-saving, which also comes with less complications and more patient acceptance.

Due to its thin design, excellent fit for narrow ridge and ensured safety due to the alveolar canal nerve bypass. Developed for single and multiple restorations.

Condensing thread



Avoiding inferior alveolar canal nerve



ROOTT C



Together with special condensing threads and embedded abutment with no microgaps, implant achieves excellent initial stability from the very beginning.

ROOTT^C



Significant time & cost saving



Immediate loadings



Excellent for narrow ridge



May avoid bone augmentation



Prosthetic variety

Cement retain with trimmable external platform, burnouts or cement-free option with patented telescopic abutments.

Telescopic



5mm
4mm
3mm



Easy management



TRS



TRS-mini

Clinical cases



By Dr. Alvaro Bastida

“FILO System is suitable in all clinical cases. Even esthetic area, narrow spaces, post-extraction and soft tissues management”



More cases





ANED



TRA



TOES



TOEA



TOE



HE



7mm



5mm



4mm



3mm

External platform



TCE
TCES
TCEXS



PCE
PCES
PCEXS

Intraoral

Scan/Impression

Extraoral



ANE



TRA



TOES



TOEA



TOE



HE



7mm



5mm



4mm



3mm


External platform




TCE
TCES
TCEXS



PCE
PCES
PCEXS



Digital workflow **ROOTT C**



Metal framework

Prosthesis with cement

Design/milling

Final structure

Telescopic abutments
External platform

Prosthesis with
telescopic solutions



+



TCE
TCES
TCEXS
PCE
PCES
PCEXS



70+

distributors



ROOTT UK

www.roottuk



<https://rootimplants.co.uk/>

Two Piece Spec



<https://rootimplants.co.uk/implant-range/>

Training

UK Training



<https://rootimplants.co.uk/training/>

Clinical cases

General



www.facebook.com/opendentalcommunity/

Facebook



www.facebook.com/rootimplants.uk/

One Piece Screw Spec



<https://rootimplants.co.uk/implant-range/>

How to use R



<https://rootimplants.co.uk/root-r-how-to-use-videos/>

Advanced



www.corticallyfixed.com/case-reports/

t: 020 3633 6150

e: info@rootimplants.co.uk

w: www.rootimplants.co.uk

l: www.linkedin.com/company/70899052



TRATE UK LTD
85 Great Portland Street
London, W1W 7LT
United Kingdom